



[Short Contents](#) | [Full Contents](#)

[Other books @ NCBI](#)

Navigation

About this book

**Appendix I.
Immunologists'
Toolbox**

Immunization.

The detection, measurement, and characterization of antibodies and their use as research and diagnostic tools.

➔ Isolation of lymphocytes.

Characterization of lymphocyte specificity, frequency, and function.

Detection of immunity *in vivo*.

Manipulation of the immune system.

Immunobiology ➔ Appendix I. Immunologists' Toolbox ➔
Isolation of lymphocytes.




Figure A.23. Peripheral blood mononuclear cells can be isolated from whole blood by Ficoll-Hypaque™ centrifugation. Diluted anticoagulated blood (left panel) is layered over Ficoll-Hypaque™ and centrifuged. Red blood cells and polymorphonuclear leukocytes or granulocytes are more dense and centrifuge through the Ficoll-Hypaque™, while mononuclear cells consisting of lymphocytes together with some monocytes band over it and can be recovered at the interface (right panel).

Search

☒ This book ☐ All books

☐ PubMed